



DHA, EPA Raw materials Analysis

1. Analytical Condition

YL 6500 Series

Oven : 170°C (for 0 min)- temp-programed 1 °C/min-225 °C(for 10 min)

Column : TR-WAX(30m*0.25mm*0.25um)

Carrier gas : He, 1.0 ml/min (Split ratio 50:1)

Injector : Capillary 270°C

Detector : FID 270°C

Sample Preconditioning : FAME

Injection Volume : 1.0ul (Liquid)

2. Chromatogram

